



User Manual

Thank you for purchasing the Digital Body Weight Scale by Etekcit®. This device can help you reach your physical goals by conveniently allowing you to measure body weight with reliable accuracy. Calculate and meet personal goals with the simple step on activation. The decorative platform accents any type of home decor with a subtle modern touch. The EB9312H makes healthy living simple with high accuracy and convenient automation.

**IMPORTANT:** Please read, understand, and comply with all the instructions provided in this manual before using the device. Failure to comply with the instructions given in this manual and/or using the device in ways other than the ones mentioned in this manual may result in serious injury and/or damage to the product.

### Safety Use & Care

- The glass platform is slippery when wet. Make sure the platform is dry before stepping on for measurement.
- Consult your primary-care doctor or physician before making changes to your diet, exercise plans and/or physical activities.
- This scale is not intended for commercial use.
- To weigh safely and accurately, please stand on the center of the platform and stand still until a result is calculated.
- **DO NOT** strike or drop the scale.
- **DO NOT** stand on the edge of the scale or jump on it.
- Always place the scale on a hard, dry and flat surface before measurement.
- Handle the scale with care, as it is a precision instrument.
- Clean the scale with a damp cloth and avoid getting water or any other liquid inside the scale.
- **DO NOT** use chemical or abrasive cleaning agents to clean the scale.
- Keep the scale in a cool environment with dry air.
- Always keep the scale in a horizontal position, do not store vertically.
- If the scale fails to turn on, check whether the battery is installed correctly. Replace the batteries if necessary.
- If an error appears on the display or if the scale does not turn off automatically, remove the batteries for 3 seconds and reinstall them.
- **DO NOT** overload the scale (Max weight 200 kg / 440 lbs.).
- If the scale is not in use for an extended period of time, remove the batteries to ensure battery conservation and to extend the scale service lifetime.

**Note:** Please be noted the scale comes with two built-in batteries and please remove the protective film before use.

### Features & Specifications

- Weight Capacity: 200 kg (440 lbs)
- Graduation: 0.1kg (0.2lbs)
- Display: LCD 4 Digits
- Battery: 2 x AAA 1.5V batteries
- Working Temperature: 10°C-35°C (50°F-95°F)
- Units: Kg / Lbs
- Auto-off: about 8 seconds
- Dimensions: 390 x 300 x 26mm
- Platform Thickness: 8mm
- Low battery indication
- Overload indication
- Re-test indication
- Backlit LCD display
- Solid glass wide platform

### Set-Up

The scale operates on the power of 2 AAA 1.5V batteries. The batteries are preinstalled in the scale. To activate, open the battery compartment cover located on the backside of the scale and remove the protective film to activate.

### Battery Replacement

1. Open the battery compartment cover on the backside of the scale.
2. Take out the used batteries and dispose them properly.
3. Install 2 new AAA 1.5V batteries, taking note of the correct polarity displayed in the battery compartment.

### Operation

1. For best accuracy, place the scale on a hard, flat floor. Avoid carpet or soft surfaces.
2. Gently step on the scale. The scale will automatically activate.
3. Remain standing evenly on the scale, without moving or shifting your feet or body while the scale processes the results. The scale will display and lock the stabilized result.

### Auto-Off

When you step off the scale, it will automatically switch off after approximately 8 seconds of inactivity.

### Converting Weight Units

This scale calculates weight measurements in both metric (kg) and Imperial (lbs) units.

- The switch to convert between units is located on the backside of the scale.
- Slide the unit switch to the desired unit at any time and the LCD screen will display the corresponding unit.

Error Message	Reason	Solution
“Err”	Scale overload. Maximum weight capacity: 200 kg, 440 lbs.	Step off the scale to avoid damage.
“Lo”	Battery power is low.	Install new batteries.
“C”	Wight reading error.	Step off the scale and re-test weight.

### Package Contents

- 1 x Digital Body Weight Scale
- 1 x User Manual
- 2 x AAA 1.5V batteries (Installed)